

Nutrimix® DC

Dual chamber bag for stable separate storage.



Clinical Nutrition

New

B | BRAUN
SHARING EXPERTISE

Nutrimix® DC

Dual chamber bag for easy and stable preparation



Filling

The Nutrimix® DC enables to compound TPN solutions with a separate storage e.g. of the lipid emulsion. This facilitates a longer shelf life of the TPN regimens if, for example, the admixtures have a limited stability due to high electrolyte concentrations or low amino acid / glucose content.

Storage

The Nutrimix® DC can be stored in the refrigerator. Before administration, the separator has to be removed. By agitating the bag, a homogeneous TPN regime is achieved.





User benefits:

- Longer storage of compounded lipid-contained TPN formulations
- Completely PVC-free bag
- Non-Latex additive port
- High flow filling port

Administration

After the mixing, the content of the Nutrimix® DC can be administered as usual.



Product specifications

Product Type	Product Specifications	Vol.	Filling lines	Units/Packages	Packages/Carton	Code No. (REF)
Nutrimix®	Nutrimix® 3/3	3000 ml	3	1	25	2112151
	Nutrimix® 2/3	2000 ml	3	1	25	2112150
	Nutrimix® 1/0	1000 ml	0	1	25	2112147
	Nutrimix® 0.5/0	500 ml	0	1	25	2112146
	Nutrimix® 0.2/0	200 ml	0	1	25	2112145
	Transfer Set 0/3	0	3	1	10	2124629
Nutrimix® (Kits of 5 bags)	Nutrimix® 3/0	3000 ml	0	5	10	2112439
	Nutrimix® 2/0	2000 ml	0	5	10	2112438
	Nutrimix® 1/0	1000 ml	0	5	10	2112437
	Nutrimix® 0.5/0	500 ml	0	5	10	2112436
	Nutrimix® 0.2/0	200 ml	0	5	10	2112435
Nutrimix® Barrier (Kits of 5 bags)	Nutrimix® Barrier 3/0	3000 ml	0	5	10	2112444
	Nutrimix® Barrier 2/0	2000 ml	0	5	10	2112443
	Nutrimix® Barrier 1/0	1000 ml	0	5	10	2112442
	Nutrimix® Barrier 0,5/0	500 ml	0	5	10	2112441
	Nutrimix® Barrier 0,2/0	200 ml	0	5	10	2112440
Nutrimix® DC	Nutrimix® dual chamber bags (2500 + 500 ml)	3000 ml	0	1	25	2112446
	Nutrimix® dual chamber bags (1250 + 250 ml)	1500 ml	0	1	25	2112445
	Transfer Set	0	4	1	25	2112449
Nutrimix® Accessories	Nutrimix® Male/Male Connector	-	-	1	50	2112447
	Nutrimix® Female/Female Connector	-	-	1	50	2112448
	Paedifix® N	0	4	1	10	2124602
	Light protection cover	-	-	-	50	2129035

**Bags 100 % PVC
and Latex free**