

Injekt 2-komponent sprøyter luer og luer-lock

Materiale:

Sprøytekammer: polypropylen

Sprøytetempel: Polyetylen



Produkt	Art.nr Luer	Art. nr: Luer Lock	Graderingsintervall
Injekt 1 ml	9166017V		0,01 ml
Injekt 2 ml	4606027V	4606701V	0,1 ml
Injekt 5 ml	4606051V	4606710V	0,2 ml
Injekt 10 ml	4606108V	4606728V	0,5 ml
Injekt 20 ml	4606205V	4606736V	1,0 ml