




Forlengesslanger

Spiralslanger	Lengde (cm)	Fyllvolum (ml)	Indre ø (mm)	Materiale	Salgsenhet stk.	Artikkel nr
	150	1,2	1,0	PUR	50	4092937
	150	4,7	2,0	PUR	50	4090365
	300	2,3	1,0	PUR	50	4092945
	300	9,4	2,0	PUR	50	4090373

Perfusor®-Slanger	Lengde (cm)	Fyllvolum (ml)	Indre ø (mm)	Materiale	Spesielle egenskaper	Salgsenhet stk.	Artikkel nr
	30	0,3	1,2	PVC	Luer-slip kobling	50	4256000
	30	0,3	1,2	PVC	Luer-Lock kobling	50	4256034
	50	0,4	1,0	PE		100	8255059
	50	0,9	1,5	PVC		100	8255172
	75	0,5	0,9	PVC		25	8722870N
	100	0,8	1,0	PE		100	8255067
	150	1,0	0,9	PVC		25	8255504N
	150	2,6	1,5	PVC		100	8722960
	150	1,2	1,0	PE		100	8722935
	150	1,2	1,0	PE	Lysbeskyttet, sort	100	8723010
	150	2,6	1,5	PVC	Lysbeskyttet, sort	100	8722919
	150	1,2	1,0	PE	Lysbeskyttet oransj transparent	100	8723017
	200	1,6	1,0	PE		100	8723060
	200	3,5	1,5	PVC		100	8722862
	200	3,8	1,5	PVC	0,2 µm filter	50	8723001
	250	2,0	1,0	PE		100	8272565
	250	4,4	1,5	PVC		100	8255490
	300	5,3	1,5	PVC		100	8255253
	300	2,4	1,0	PE		100	8250146

Heidelberger forlengesslanger	Lengde (cm)	Fyllvolum (ml)	Indre ø (mm)	Materiale	Salgsenhet stk.	Artikkel nr
	30	2,1	3,0	PVC	100	4097300
	75	5,3	3,0	PVC	100	4097173
	100	7,1	3,0	PVC	100	4097262
	150	10,6	3,0	PVC	100	4097408
	200	14,1	3,0	PVC	100	4097190
	320	22,6	3,0	PVC	100	4097130