





# MOBILITY AND SAFETY IN A NEW DIMENSION



COLOR-CODED SPIRAL LINES BY B. BRAUN

Single lumen	Inner Ø	Lengths (cm)	Color	Item number	Pcs./car
	1 mm	150	<input type="checkbox"/>	4092937	50
	1 mm	300	<input type="checkbox"/>	4092945	50
	1 mm	500	<input type="checkbox"/>	4092953	50
	1 mm	700	<input type="checkbox"/>	4092961	50
	1 mm	850	<input type="checkbox"/>	4092970	50
	2 mm	150	<input type="checkbox"/>	4090365	50
	2 mm	300	<input type="checkbox"/>	4090373	50
	2 mm	500	<input type="checkbox"/>	4090438	50
	2 mm	700	<input type="checkbox"/>	4090381	50
	2 mm	850	<input type="checkbox"/>	4090390	50

Double lumen	Inner Ø	Lengths (cm)	Color	Item number	Pcs./car
	1 mm	300	<input checked="" type="checkbox"/> <input type="checkbox"/>	4091736	50
	1 mm	300	<input type="checkbox"/> <input checked="" type="checkbox"/>	4091752	50
	1 mm	700	<input type="checkbox"/> <input checked="" type="checkbox"/>	4093607	50
	1 mm	850	<input type="checkbox"/> <input checked="" type="checkbox"/>	4091728	50
	2 mm	300	<input checked="" type="checkbox"/> <input type="checkbox"/>	4093830	50
	2 mm	500	<input checked="" type="checkbox"/> <input type="checkbox"/>	4093850	40
	2 mm	700	<input type="checkbox"/> <input checked="" type="checkbox"/>	4093870	30
	2 mm	850	<input type="checkbox"/> <input checked="" type="checkbox"/>	4093885	25

Triple lumen	Inner Ø	Lengths (cm)	Color	Item number	Pcs./car
	1 mm	150	■ ■ ■	4093115	50
	1 mm	300	■ ■ ■	4093130	50
	1 mm	500	■ ■ ■	4093150	50
	1 mm	700	■ ■ ■	4093170	50
	1 mm	850	■ ■ ■	4093185	50
	2 mm	150	■ ■ ■	4093215	50
	2 mm	300	■ ■ ■	4093230	50
	2 mm	500	■ ■ ■	4093250	40
	2 mm	700	■ ■ ■	4093270	25
	2 mm	850	■ ■ ■	4093285	25

Lengths	150 cm	300 cm	500 cm	700 cm	850 cm
Diameter 2 x 3.2 mm/Filling Volume	4.7 ml	9.9 ml	15.7 ml	22 ml	26.7 ml
Flow Rate ml/min (BS)	ca. 61	ca. 40	ca. 28	ca. 20	ca. 17
Diameter 1 x 1.9 mm/Filling Volume	1.2 ml	2.3 ml	3.9 ml	5.5 ml	6.6 ml
Flow Rate ml/min (BS)	ca. 6.0	ca. 3.5	ca. 2.5	ca. 1.7	ca. 1.5

#### PROSET – THE SET FOR COMPLEX INFUSION THERAPY

Our sales representatives are available to advise you regarding the choice of the most suitable equipment for your specific processes and clinical requirements.

The new color-coded spiral lines are the central element for all functions, including the Discofix® C Manifolds distribution system, the filtering and the patient access.

